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
Isocratic Separation of Anilines on HALO C18

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

TEST CONDITIONS:

Column: HALO C18, 4.6 x 50 mm
 Part Number: 92814-402
 Mobile Phase: 60/40– A/B
 A=0.02 M sodium phosphate buffer, pH=7.0
 B= Acetonitrile
 Flow Rate: 2.0 mL/min.
 Pressure: 211 Bar
 Temperature: 25°C
 Detection: UV 254 nm, VWD
 Injection Volume: 1.0 µL
 Sample Solvent: ACN/ water-50/50
 Response Time: 0.02 sec.
 Flow Cell: 2.5 µL semi-micro
 LC System: Shimadzu Prominence UFLC XR
 Extra column volume: ~14 µL


STRUCTURES:


p-Aminobenzoic acid


3-Nitroaniline


1,2-phenylenediamine


2-Nitroaniline

p-Anisidine


4-Chloroaniline

Aniline

Application Note: 35-EX

Isocratic Separation of Dinitrotoluenes on HALO RP-Amide Phase HALO RP-Amide 色谱柱快速等度分离二硝基甲苯类化合物


PEAK IDENTITIES:

1. Uracil (尿嘧啶)
2. 2,4-dinitrotoluene (2,4-二硝基甲苯)
3. 2,6-dinitrotoluene (2,6-二硝基甲苯)
4. 3,4-dinitrotoluene (3,4-二硝基甲苯)
5. 2,3-dinitrotoluene (2,3-二硝基甲苯)


TEST CONDITIONS:

Column: 4.6 x 50 mm, HALO RP-Amide
Part Number: 92814-407
Mobile Phase: 80/20-Water/Acetonitrile
Flow Rate: 2.5 mL/min.
Pressure: 257 Bar
Temperature: 27 °C
Detection: UV 254 nm, VWD
Injection Volume: 1.0 µL
Sample Solvent: 50/50-Methanol/Acetonitrile
Response Time: 0.02 sec.
Flow Cell: 2.5 µL semi-micro
LC System: Shimadzu Prominence UFLC XR
Extra column volume: ~14 µL


STRUCTURES


Uracil


2,6-Dinitrotoluene

2,4-Dinitrotoluene

3,4-Dinitrotoluene


2,3-Dinitrotoluene

These dinitrotoluene isomers are difficult to separate, but can be separated with near baseline resolution in under 7 minutes using a HALO RP-Amide column.

Application Note: 53-G

Separation of Mixed Polarity Compounds on HALO C18 and ES-CN

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






1-Chloro-4-nitrobenzene




Toluene

These separations of polar and non-polar compounds show significant differences in selectivity between HALO C18 and ES-CN stationary phases. Note the increased retention of nitro compounds and reduced retention of non-polar compounds on HALO ES-CN phase.

Application Note: 23-N

Separation of Neutral Aromatics on HALO PFP, C18 and Phenyl-Hexyl HALO PFP, C18 和苯基色谱柱分离芳香化合物对比


PEAK IDENTITIES:

1. Butylbenzene (丁基苯)
 2. Acenaphthene (萵)
 3. 1-Phenylnaphthalene (1-苯基萘)
 4. Pyrene (苊)
- i=impurities


TEST CONDITIONS:

Column: 4.6 x 50 mm, HALO PFP, C18, Phenyl-Hexyl
 Part Numbers: 92814-409, -402, -406, resp.
 Mobile Phase: 30/70-water/methanol
 Flow Rate: 2.0 mL/min.
 Pressure: approximately 250 Bar
 Temperature: 40 °C
 Detection: UV 254 nm, VWD
 Injection Volume: 1.0 µL
 Sample Solvent: methanol
 Response Time: 0.02 sec.
 Flow Cell: 2.5 µL semi-micro
 LC System: Shimadzu Prominence UFLC XR
 Extra column volume: ~14 µL


STRUCTURES:

Butylbenzene

1-Phenylnaphthalene

Acenaphthene


Pyrene

The separation of non-polar aromatic compounds on these three Halo bonded phases under the same conditions show differences in selectivity that can be utilized in optimizing difficult separations.

Application Note: 085-PS

Separation of Galanthamine and Quetiapine on HALO PFP (pentafluorophenyl) Phase

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

Galanthamine


Quetiapine

These psychiatric drugs can be rapidly separated on a HALO PFP column in just one minute.

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

Hesperidin

Diosmin


These two semisynthetic flavonoids can be rapidly separated using HALO-5 PFP (pentafluorophenyl) stationary phase at a low pressure. Note that just the addition of a double bond results in a difference that allows these two very similar compounds to be separated.

*Needed for solubility reasons.

Application Note: 083-FL

Separation of Diosmin and Hesperidin on HALO Phenyl-Hexyl and HALO RP-Amide

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
1. Diosmin (香叶木甙)
2. Hesperidin (橙皮素)


TEST CONDITIONS:

Column 1: 4.6 x 50 mm, HALO Phenyl-Hexyl
 Part Number: 92814-406
 Column 2: 4.6 x 50 mm, HALO RP Amide
 Part Number: 92814-407
 Mobile Phase: 78/22: Water/acetonitrile
 Flow Rate: 1.5 mL/min.
 Pressure: 145 Bar
 Temperature: 40°C
 Detection: UV 254 nm, VWD
 Injection Volume: 0.5 µL
 Sample Solvent: Dimethylformamide*
 Response Time: 0.02 sec.
 Flow Cell: 2.5 µL semi-micro
 LC System: Shimadzu Prominence UFLC XR
 ECV: ~14 µL

STRUCTURES:

Diosmin

Hesperidin


These two semi-synthetic flavonoids are often taken to enhance vascular health. The two compounds may be easily separated using either HALO RP-Amide or HALO Phenyl-Hexyl phases. Note the difference in elution order on the two phases.

*Needed for solubility reasons.

Application Note: 24-P

Isocratic Separation of Phthalate Esters on HALO C18

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

TEST CONDITIONS:

Column: 4.6 x 50 mm, HALO C18
 Part Number: 92814-402
 Mobile Phase: 20/80 A/B
 A=water
 B= acetonitrile
 Flow Rate: 1.5 mL/min.
 Pressure: 97 Bar
 Temperature: 27°C
 Detection: UV 254 nm, VWD
 Injection Volume: 0.5 µL
 Sample Solvent: acetonitrile
 Response Time: 0.02 sec.
 Flow Cell: 2.5 µL semi-micro
 LC System: Shimadzu Prominence UFLC XR
 Extra column volume: ~14 µL


STRUCTURES:


Uracyl


Di-n-propylphthalate

Dimethylphthalate

Di-n-butylphthalate


Diethylphthalate

Application Note: 61-CB

Separation of Carbamate Pesticides on HALO C18 Phase

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


Prothion


Chlorprothion

This separation illustrates a rapid HPLC determination of three carbamate pesticides on HALO C18 phase in just under a minute. The unique Fused-Core technology allows the use of high flow rates at moderate pressures while retaining high efficiency.

Application Note: 55-PU

Separation of Phenyl Urea Pesticides on HALO Phenyl-Hexyl Phase HALO 苯基 - 己基相快速分离苯基脲农药


PEAK IDENTITIES:

1. Fenuron (非草隆)
2. Monuron (灭草隆)
3. Fluomethuron (伏草隆)
4. Isoproturon (异丙隆)
5. Diuron (敌草隆)
6. Siduron A (环草隆A)
7. Siduron B (环草隆B)
8. Linuron (利谷隆)
9. Neburon (草不隆)


TEST CONDITIONS:

Column: 4.6 x 100 mm, HALO Phenyl-Hexyl
Part Number: 92814-606
Mobile Phase: 50/50-A/B
A= 0.025 M Potassium phosphate buffer, adj. to pH = 2.5
B=Acetonitrile
Flow Rate: 1.5 mL/min.
Pressure:220 Bar
Temperature: 30 °C
Detection: UV 245 nm, VWD
Injection Volume: 0.5 µL
Sample Solvent: Acetonitrile
Response Time: 0.02 sec.
Flow Cell: 2.5 µL semi-micro
LC System: Shimadzu Prominence UFLC XR


STRUCTURES:


Fenuron


Isoproturon


Siduron B


Monuron


Diuron


Linuron

Fluomethuron

Siduron A


Neburon

This separation illustrates the use of the highly efficient HALO Fused-Core Phenyl-Hexyl stationary phase in the analysis of common herbicides. The short run times allow analyses using isocratic conditions so that column equilibration time is not required between runs.

Application Note: 59-PU

Separation of Phenyl Urea Pesticides on HALO C18 Phase

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

TEST CONDITIONS:

Column: 4.6 x 100 mm, HALO C18
 Part Number: 92814-602
 Mobile Phase: 50/50-A/B
 A= 0.025 M Potassium phosphate buffer, adj. to pH = 2.5
 B=Acetonitrile
 Flow Rate: 2.0 mL/min.
 Pressure: 300 Bar
 Temperature: 30 °C
 Detection: UV 245 nm, VWD
 Injection Volume: 0.5 µL
 Sample Solvent: Acetonitrile
 Response Time: 0.02 sec.
 Flow Cell: 2.5 µL semi-micro
 LC System: Shimadzu Prominence UFLC XR


STRUCTURES:


Fenuron


Isoproturon


Siduron B


Monuron


Diuron


Linuron

Fluomethuron

Siduron A


Neburon

This separation illustrates the use of the highly efficient 2.7 µm HALO Fused-Core C18 stationary phase in the analysis of some popular herbicides. The short run times allow analyses using isocratic conditions so that column equilibration time is not required between runs.

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




1. Amoxicillin (阿莫西林)
2. Ampicillin (氨苄西林)


TEST CONDITIONS:

Column: 4.6 x 50 mm, HALO RP-Amide
Part Number: 92814-407
Mobile Phase: 82/18: A/B
A= 0.02 M Phosphate buffer, pH=2.7
B= Acetonitrile
Flow Rate: 2.0 mL/min.
Pressure: 200 Bar
Temperature: 30°C
Detection: UV 240 nm, VWD
Injection Volume: 1.0 µL
Sample Solvent: 80/20: water/acetonitrile
Response Time: 0.02 sec.
Flow Cell: 2.5 µL semi-micro
LC System: Shimadzu Prominence UFLC XR
ECV: ~14 µL

STRUCTURES:

Amoxicillin


Ampicillin

This chromatogram shows a time-efficient separation of two weakly retained antibiotics on a short HALO RP-Amide column.

Application Note: 58-AM

Isocratic Separation of Amphenicols on HALO RP-Amide Phase

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

Thiamphenicol


Chloramphenicol

This separation shows a rapid HPLC method for the analysis of amphenicols using HALO RP-Amide phase. To improve the sensitivity of detection the first peak was monitored @ 240 nm and the second @ 280 nm.

Application Note: 080-AF

Separation of Antibiotic and Antifungal Drugs on HALO RP-Amide

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Detection: UV 230 nm, VWD

Injection Volume: 0.3 µL


Sample Solvent: water/acetonitrile-25/75

Response Time: 0.02 sec.


Flow Cell: 2.5 µL semi-micro

LC System: Shimadzu Prominence UFLC XR


STRUCTURES:


Ketoconazole


Econazole


Naftifine


Sulconazole

Clotrimazole

Tolnaftate


Clofazimine

The antimicrobial drug clofazimine and the antifungal drugs can be rapidly analyzed using the HALO RP-Amide column under gradient conditions with low back pressure.

Application Note: 067-AD

Separation of Antidepressants on HALO Penta-HILIC HALO Penta-HILIC 色谱柱快速分离抗抑郁类药物

PEAK IDENTITIES:


1. Trimipramine (三甲丙咪嗪)
2. Amitriptyline (阿米替林)
3. Doxepin (多虑平)
4. Nortriptyline (去甲替林)
5. Amoxapine (阿莫沙平)

TEST CONDITIONS:


Column: 4.6 x 100 mm, HALO Penta-HILIC
 Part Number: 92814-605
 Mobile Phase: 7/93: A/B
 A= 0.1 M Ammonium formate, pH=3.5 (adj.)
 B= Acetonitrile
 Flow Rate: 2.5 mL/min.
 Pressure: 165 Bar
 Temperature: 30°C
 Detection: UV 254 nm, VWD
 Injection Volume: 0.5 µL
 Sample Solvent: 10/90: Water/Acetonitrile
 Response Time: 0.02 sec.
 Flow Cell: 2.5 µL semi-micro
 LC System: Shimadzu Prominence UFLC XR
 ECV: ~14 µL

Basic drugs such as antidepressants can be separated rapidly under HILIC conditions with good peak shape using HALO Penta-HILIC stationary phase.


STRUCTURES:


Trimipramine


Doxepin

Amitriptyline

Nortriptyline


Amoxapine

HPLC Application Note: 72-AB

Separation of Penicillins on HALO Phenyl-Hexyl

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

TEST CONDITIONS:

Column: 4.6 x 50 mm, HALO Phenyl-Hexyl
 Part Number: 92814-406
 Mobile Phase: 40/60: A/B
 A= 0.02 M Phosphate buffer, pH=3.0 (adj.)
 B= Methanol
 Flow Rate: 1.5 mL/min.
 Pressure: 200 Bar
 Temperature: 40°C
 Detection: UV 230 nm, VWD
 Injection Volume: 1.0 µL
 Sample Solvent: 50/50 Water/acetonitrile
 Response Time: 0.02 sec.
 Flow Cell: 2.5 µL semi-micro
 LC System: Shimadzu Prominence UFLC XR
 ECV: about 14 µL


STRUCTURES:


Penicillin G

Oxacillin

Piperacillin


Cloxacillin

These four penicillin drugs can be rapidly separated on HALO Fused-Core Phenyl-Hexyl bonded phase columns.

HPLC Application Note: 69-AB

Separation of Cephalosporins on HALO ES-CN HALO ES-CN 色谱柱快速分离头孢菌素类药物


PEAK IDENTITIES:

1. Cefadroxil (头孢羟氨苄)
2. Ceftazidime (头孢他啶)
3. Cefaclor (头孢克洛)
4. Cephalexin (头孢氨苄)
5. Cephadrine (头孢拉定)
6. Cefotaxime (头孢噻肟)
7. Cefoxitin (头孢西丁)
8. Cefazolin (头孢唑啉)
9. Cephalothin (头孢菌素)


TEST CONDITIONS:

Columns: HALO ES-CN, 4.6 x 50 mm
Part Numbers 92814-404
Mobile Phase: A= 0.02 M Phosphate buffer, pH=2.7 (adj.)
B= Methanol
Gradient: 20% B to 40% B in 2.5 min.
Flow Rate: 2.0 mL/min.
Pressure: 225 bar at start of gradient
Temperature: 40° C
Detection: UV 254 nm, VWD
Injection Volume: 1.0 µL
Sample Solvent: 70/30 Water/Methanol
Response Time: 0.02 sec.
Flow Cell: 2.5 µL semi-micro
LC System: Shimadzu Prominence UFLC XR
ECV: about 14 µL


STRUCTURES:


Cefadroxil


Cephalexin


Cefoxitin


Ceftazidime


Cephadrine


Cefazolin

Cefaclor

Cefotaxime


Cephalothin

These cephalosporins can be rapidly separated by reversed-phase HPLC on a HALO Fused-Core ES-CN bonded phase column.

Application Note: 078-DU

HPLC Separation of Diuretics on HALO Phenyl-Hexyl HALO 苯基色谱柱快速分离利尿剂类药物

PEAK IDENTITIES:

1. Amiloride (阿米洛利)
2. Caffeine (咖啡因)
3. Chlorothiazide (利噻嗪)
4. Hydrochlorothiazide (氢氯噻嗪)
5. Triamterene (氨苯蝶啶)
6. Torsemide (托拉塞米)
7. Furosemide (呋喃苯胺酸)
8. Indapamide (吲达帕胺)
9. Bumetanide (布美他尼)

TEST CONDITIONS:

Column: 4.6 x 100 mm, HALO Phenyl-Hexyl
Part Number: 92814-606
Mobile Phase:
A= 0.02 M Potassium phosphate buffer, (pH=3)
B= Acetonitrile

Time	%B	Time	%B
0	15	3.0	50
1.7	15	7.0	60

Flow Rate: 1.5 mL/min.

Pressure: 253 Bar (@start of gradient)

Temperature: 30°C

Detection: UV 230 nm, VWD

Injection Volume: 2 µL

Sample Solvent: Acetonitrile


Response Time: 0.02 sec.

Flow Cell: 2.5 µL semi-micro


LC System: Shimadzu Prominence UFLC XR

ECV: ~14 µL


STRUCTURES:


Amiloride


Hydrochlorothiazide


Furosemide


Caffeine


Triamterene


Indapamide

Chlorothiazide

Torsemide


Bumetanide

This separation illustrates the utility of HALO Fused-Core Phenyl-Hexyl phase in the rapid analysis of common diuretics.

Application Note: 14-NS

Gradient Separation of NSAIDS on HALO C8 HALO C8 色谱柱快速梯度分离非甾体类药物

PEAK IDENTITIES:

1. Acetaminophen (对乙酰氨基酚)
2. Aspirin (阿司匹林)
3. Salicylic acid (水杨酸)
4. Tolmetin (甲苯酰吡啶乙酸)
5. Ketoprofen (酪洛芬)
6. Naproxen (萘普生)
7. Fenoprofen (非诺洛芬)
8. Diclofenac (双氯芬酸)
9. Ibuprofen (布洛芬)
10. Mefenamic acid (甲灭酸)

TEST CONDITIONS:

Column: 4.6 x 50 mm, HALO C8

Part Number: 92814-408

Mobile Phase: 38/62-- A/B (start)

A= 0.02 M sodium phosphate buffer, pH=2.5

B= methanol

Gradient:

0.0 min 62%B

0.1 min 62%B

2.0 min 85%B

Flow Rate: 2.0 mL/min.

Pressure: 286 Bar

Temperature: 35°C

Detection: UV 254 nm, VWD

Injection Volume: 1.0 µL

Sample Solvent: mobile phase


Response Time: 0.02 sec.

Flow Cell: 2.5 µL semi-micro


LC System: Shimadzu Prominence UFLC XR

Extra column volume: ~14 µL


STRUCTURES:


Acetaminophen


Aspirin


Salicylic acid


Tolmetin


Ketoprofen


Naproxen


Fenoprofen

Diclofenac

Ibuprofen


Mefenamic acid

Application Note: 56-NS

Isocratic Separation of NSAIDs on HALO ES-CN Phase

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








Fenoprofen


Mefenamic acid

This separation illustrates the separating power of HALO Fused-Core stationary phases. Nine NSAID drugs are separated in under one minute on a 50 mm HALO ES-CN column.

Application Note: 16-NS

Gradient Separation of NSAIDS on HALO RP-Amide HALO RP-Amide 色谱柱梯度分离非甾体抗炎药

PEAK IDENTITIES:

1. Acetaminophen (对乙酰氨基酚)
2. Aspirin (阿司匹林)
3. Salicylic acid (水杨酸)
- i=impurity
4. Tolmetin (甲苯吡啶乙酸)
5. Ketoprofen (萘洛芬)
6. Naproxen (萘普生)
7. Fenoprofen (非诺洛芬)
8. Ibuprofen (布洛芬)
9. Diclofenac (双酚氯酸)
10. Mefenamic acid (甲灭酸)

TEST CONDITIONS:

Column: 4.6 x 50 mm, HALO RP-Amide
Part Number: 92814-407
Mobile Phase: 50/50– A/B (start)
A= 0.02 M sodium phosphate buffer, pH=2.5
B= methanol

Gradient:

0.0 min 50%B
0.1 min 50%B
0.5 min 55%B
3.5 min 80%B
4.0 min 80%B

Flow Rate: 2.0 mL/min.

Pressure: 289 Bar

Temperature: 35°C

Detection: UV 254 nm, VWD

Injection Volume: 1.0 µL

Sample Solvent: mobile phase


Response Time: 0.02 sec.

Flow Cell: 2.5 µL semi-micro


LC System: Shimadzu Prominence UFLC XR

Extra column volume: ~14 µL


STRUCTURES:


Acetaminophen


Naproxen


Aspirin


Fenoprofen


Salicylic acid


Diclofenac


Tolmetin

Ibuprofen



Ketoprofen

Mefenamic acid

HALO-5® C18 Lot to Lot Reproducibility HALO-5 C18 色谱柱批次重复性实验结果

PEAK IDENTITIES:

1. Uracil (尿嘧啶)
2. Phenol (苯酚)
3. 4-Cl-1-Nitrobenzene (4-氯硝基苯)
4. Naphthalene (萘)

TEST CONDITIONS:

Columns: 4.6 x 50 mm HALO-5 C18
Part Number: 95814-402
Mobile Phase: 57/43 acetonitrile/water
Flow Rate: 1.0 mL/min.
Pressure: 39 bar
Temperature: 30° C
Detection: UV 254 nm, VWD
Injection Volume: 1.0 µL
Sample Solvent: 50/50 ACN/water
Flow Cell: 5 µL semi-micro
LC System: Agilent 1100

The retention factor and selectivity calculated across several batches of HALO-5 C18 show superior reproducibility. Retention factor is calculated for naphthalene while selectivity is calculated between naphthalene and 4-chloro-1-nitrobenzene.

References

Fused-Core particle technology has been thoroughly studied by recognized experts in the field of chromatography. Below is a partial list of their publications that report on the characteristics, performance, and application of UHPLC columns packed with Fused-Core particles.

1. Mass-transfer kinetics in a shell packing material for chromatography. *Anal. Chem* (2007) 79, 5972-5979, Cavazzini, A., Gritti, F., Kaczmarek, K., Marchetti, N. & Guiochon, G.
2. Fused core particle technology as an alternative to sub 2 m particles to achieve high separation efficiency with low backpressure. *Journal of separation science* (2007) 30, 3104-3109, Cunliffe, J. & Maloney, T.
3. Comparison between the efficiencies of columns packed with fully and partially porous C18-bonded silica materials. *Journal of Chromatography A* (2007) 1157, 289-303, Gritti, F., Cavazzini, A., Marchetti, N. & Guiochon, G.
4. Unusual behavior of the height equivalent to a theoretical plate of a new poroshell stationary phase at high temperatures. *Journal of Chromatography A* (2007) 1169, 125-138, Gritti, F. & Guiochon, G.
5. Comparative study of the performance of columns packed with several new fine silica particles: Would the external roughness of the particles affect column properties? *Journal of Chromatography A* (2007) 1166, 30-46, Gritti, F. & Guiochon, G.
6. Comparison between the loading capacities of columns packed with partially and totally porous fine particles: What is the effective surface area available for adsorption? *Journal of Chromatography A* (2007) 1176, 107-122, Gritti, F. & Guiochon, G.
7. Fused-core silica column high-performance liquid chromatography/tandem mass spectrometric determination of rimonabant in mouse plasma. *Anal. Chem* (2007) 79, 5668-5673, Hsieh, Y., Duncan, C. & Brisson, J.
8. Modeling of the mass-transfer kinetics in chromatographic columns packed with shell and pellicular particles. *Anal. Chem* (2007) 79, 4648-4656, Kaczmarek, K. & Guiochon, G.
9. Fused core particles for HPLC columns. *American laboratory* (2007) 39, 18-21, Kirkland, J., Langlois, T. & DeStefano, J.
10. New Chromatography Columns and Accessories at Pittcon 2007: Part I. LC GC North America (2007) March, 248-266, Majors, R. E.
11. Gradient elution separation and peak capacity of columns packed with porous shell particles. *Journal of Chromatography A* (2007) 1163, 203-211, Marchetti, N., Cavazzini, A., Gritti, F. & Guiochon, G.
12. High peak capacity separations of peptides in reversed-phase gradient elution liquid chromatography on columns packed with porous shell particles. *Journal of Chromatography A* (2007) 1176, 206-216, Marchetti, N. & Guiochon, G.
13. Kinetic plot method as a tool to design coupled column systems producing 100,000 theoretical plates in the shortest possible time. *Journal of Chromatography A* (2008) 1212, 23-34, Cabooter, D., Lestremieu, F., Lynen, F., Sandra, P. & Desmet, G.
14. Determination of zearalenone and its metabolites in urine samples by liquid chromatography with electrochemical detection using a carbon nanotube-modified electrode. *Journal of Chromatography A* (2008) 1212, 54-60, De Andres, F., Zougagh, M., Casta eda, G. & Ríos, A.
15. Characteristics of Superficially-Porous Silica Particles for Fast HPLC: Some Performance Comparisons with Sub-2-µm Particles. *Journal of chromatographic science* (2008) 46, 254-260, DeStefano, J., Langlois, T. & Kirkland, J.
16. Surface Area of Reversed-Phase HPLC Columns. *Analytical chemistry* (2008) 80, 6358-6364, Giaquinto, A., Liu, Z., Bach, A. & Kazakevich, Y.
17. Implementations of two-dimensional liquid chromatography. *Journal of Chromatography A* (2008) 1189, 109-168, Guiochon, G., Marchetti, N., Mriziq, K. & Shalliker, R.
18. Serial coupled columns reversed-phase separations in high-performance liquid chromatography: Tool for analysis of complex real samples. *Journal of Chromatography A* (2008) 1188, 208-215, Herrero, M., Cacciola, F., Donato, P., Giuffrida, D., Dugo, G., Dugo, P. & Mondello, L.
19. RPLC-ion-trap-FTMS method for lipid profiling of plasma: method validation and application to p53 mutant mouse model. *Journal of Proteome Research* (2008) 7, 4982-4991, Hu, C., van Dommelen, J., van der Heijden, R., Spijksma, G., Reijmers, T., Wang, M., Slee, E., Lu, X., Xu, G. & van der Greef, J.
20. High-performance liquid chromatography-atmospheric pressure photoionization/tandem mass spectrometry for the detection of 17α-ethinylestradiol in hepatocytes. *Journal of Chromatography B* (2008) 870, 186-191, Li, F., Hsieh, Y. & Korfmacher, W.
21. Evaluation of the properties of a superficially porous silica stationary phase in hydrophilic interaction chromatography. *Journal of Chromatography A* (2008) 1193, 85-91, McCalley, D.
22. Estimation of the extent of the water-rich layer associated with the silica surface in hydrophilic interaction chromatography. *Journal of Chromatography A* (2008) 1192, 225-229, McCalley, D. & Neue, U.
23. Determination of fatty alcohol ethoxylates by derivatization with phthalic anhydride followed by liquid chromatography with UV-vis detection. *Journal of Chromatography A* (2008) 1203, 47-53, Micó-Tormos, A., Simó-Alfonso, E. & Ramis-Ramos, G.
24. Synthesis, Structure- Activity Relationships, and Biological Profiles of a Quinazolinone Class of Histamine H3 Receptor Inverse Agonists. *Journal of medicinal chemistry* (2008) 51, 4780-4789, Nagase, T., Mizutani, T., Ishikawa, S., Sekino, E., Sasaki, T., Fujimura, T., Ito, S., Mitobe, Y., Miyamoto, Y. & Yoshimoto, R.
25. Synthesis and evaluation of structurally constrained quinazolinone derivatives as potent and selective histamine H3 receptor inverse agonists. *J. Med. Chem* (2008) 51, 6889-6901, Nagase, T., Mizutani, T., Sekino, E., Ishikawa, S., Ito, S., Mitobe, Y., Miyamoto, Y., Yoshimoto, R., Tanaka, T. & Ishihara, A.
26. Fused-core particles: a practical alternative to sub-2 micron particles. *Journal of chromatographic science* (2008) 46, 883-886, Salisbury, J.
27. What is "Hot" in Column Technologies for Liquid Chromatography? *American laboratory* (2008) 40, 13-17, Svec, F.
28. Particle packed columns and monolithic columns in high-performance liquid chromatography-comparison and critical appraisal. *Journal of Chromatography A* (2008) 1184, 393-415, Unger, K., Skudas, R. & Schulte, M.

References

29. Radial heterogeneity of some analytical columns used in high-performance liquid chromatography. *Journal of Chromatography A* (2009) 1216, 3185-3191, Abia, J., Mriziq, K. & Guiochon, G.
30. Application of a Dried-DMSO rapid throughput 24-h equilibrium solubility in advancing discovery candidates. *European Journal of Pharmaceutical Sciences* (2009) 37, 172-182, Alelyunas, Y., Liu, R., Pelosi-Kilby, L. & Shen, C.
31. Comprehensive two-dimensional liquid chromatography separations of pharmaceutical samples using dual Fused-Core columns in the 2nd dimension. *Journal of Chromatography A* (2009) 1216, 1338-1345, Alexander, A. & Ma, L.
32. Looking Beyond the Column: An Investigation into Method Irreproducibility for an Active Pharmaceutical Ingredient. *Chromatographia* (2009) 70, 1561-1567, Bernardoni, F., Sajonz, P., Zang, J., Lee, C., Marcinko, S., Abraham, A. & Helmy, R.
33. Fused core, sub 2 µm packings, and monolithic HPLC columns: a comparative evaluation. *Journal of separation science* (2009) 32, 2723-2731, Brice, R., Zhang, X. & Colón, L.
34. A high-throughput LC-MS/MS method for the quantitation of posaconazole in human plasma: Implementing fused core silica liquid chromatography. *Journal of pharmaceutical and biomedical analysis* (2009) 50, 46-52, Cunliffe, J., Noren, C., Hayes, R., Clement, R. & Shen, J.
35. Approaches to comprehensive multidimensional liquid chromatography systems. *Journal of Chromatography A* (2009) 1216, 1363-1371, Fairchild, J., Horváth, K. & Guiochon, G.
36. Tryptic digest analysis by comprehensive reversed phase two reversed phase liquid chromatography (RP LCX 2RP LC) at different pH's. *Journal of separation science* (2009) 32, 1137-1144, François, I., Cabooter, D., Sandra, K., Lynen, F., Desmet, G. & Sandra, P.
37. Semi automated tandem mass spectrometric (MS/MS) triple quadrupole operating parameter optimization for high throughput MS/MS detection workflows. *Rapid Communications in Mass Spectrometry* (2009) 23, 1303-1312, Geddes, K., Adamson, G., Dube, N., Crathern, S. & King, R.
38. Comparison of molecularly imprinted, mixed-mode and hydrophilic balance sorbents performance in the solid-phase extraction of amphetamine drugs from wastewater samples for liquid chromatography-tandem mass spectrometry determination. *Journal of Chromatography A* (2009) 1216, 8435-8441, González-Mariño, I., Quintana, J., Rodríguez, I., Rodil, R., González-Peñas, J. & Cela, R.
39. Peak shapes of acids and bases under overloaded conditions in reversed-phase liquid chromatography, with weakly buffered mobile phases of various pH: A thermodynamic interpretation. *Journal of Chromatography A* (2009) 1216, 63-78, Gritti, F. & Guiochon, G.
40. Adsorption mechanism of acids and bases in reversed-phase liquid chromatography in weak buffered mobile phases designed for liquid chromatography/mass spectrometry. *Journal of Chromatography A* (2009) 1216, 1776-1788, Gritti, F. & Guiochon, G.
41. Evaluation of The Peak Capacity of Various RP-Columns for Small Molecule Compounds in Gradient Elution. *Chromatographia* (2009) 70, 1045-1054, Guo, Y., Srinivasan, S. & Gaiki, S.
42. Optimization strategies for off-line two-dimensional liquid chromatography. *Journal of Chromatography A* (2009) 1216, 2511-2518, Horváth, K., Fairchild, J. & Guiochon, G.
43. Hydrophilic interaction liquid chromatography/tandem mass spectrometry for the simultaneous determination of dasatinib, imatinib and nilotinib in mouse plasma. *Rapid Communications in Mass Spectrometry* (2009) 23, 1364-1370, Hsieh, Y., Galviz, G., Zhou, Q. & Duncan, C.
44. Analysis of impurities in streptomycin and dihydrostreptomycin by hydrophilic interaction chromatography/electrospray ionization quadrupole ion trap/time of flight mass spectrometry. *Rapid Communications in Mass Spectrometry* (2009) 23, 907-914, Kawano, S.
45. Characterization of the alcoholic fraction of vegetable oils by derivatization with diphenic anhydride followed by high-performance liquid chromatography with spectrophotometric and mass spectrometric detection. *Journal of Chromatography A* (2009) 1216, 230-236, Lerma-García, M., Ramis-Ramos, G., Herrero-Martínez, J., Gimeno-Adelantado, J. & Simó-Alfonso, E.
46. Equation for Peak Capacity Estimation in Two-Dimensional Liquid Chromatography. *Analytical chemistry* (2009) 81, 845-850, Li, X., Stoll, D. & Carr, P.
47. Chromatographic and spectroscopic studies on the chiral recognition of sulfated [beta]-cyclodextrin as chiral mobile phase additive:: Enantiomeric separation of a chiral amine. *Journal of Chromatography A* (2009) 1216, 1232-1240, Ma, S., Shen, S., Haddad, N., Tang, W., Wang, J., Lee, H., Yee, N., Senanayake, C. & Grinberg, N.
48. Determination of fatty alcohol ethoxylates with diphenic anhydride derivatization and liquid chromatography with spectrophotometric detection:: A comparative study with other anhydrides. *Journal of Chromatography A* (2009) 1216, 3023-3030, Micó-Tormos, A., Bianchi, F., Simó-Alfonso, E. & Ramis-Ramos, G.
49. A review of current trends and advances in modern bio-analytical methods: Chromatography and sample preparation. *Analytica chimica acta* (2009) 656, 8-35, Nováková, L. & Vlcková, H.
50. Comparison of fused-core and conventional particle size columns by LC-MS/MS and UV: Application to pharmacokinetic study. *Journal of pharmaceutical and biomedical analysis* (2009) 50, 491-500, Song, W., Pabbisetty, D., Groeber, E., Steenwyk, R. & Fast, D.
51. Reductive Amination of Aldehyde 2, 4-Dinitrophenylhydrazones Using 2-Picoline Borane and High-Performance Liquid Chromatographic Analysis. *Analytical chemistry* (2009) 81, 485-489, Uchiyama, S., Inaba, Y., Matsumoto, M. & Suzuki, G.
52. Analysis of human antibody IgG2 domains by reversed-phase liquid chromatography and mass spectrometry. *Journal of Chromatography B* (2009) 877, 1613-1620, Yan, B., Eris, T., Yates, Z., Hong, R. W., Steen, S., Kleemann, G., Wang, W. & Liu, J. L.
53. Critical comparison of performances of superficially porous particles and sub-2 µm particles under optimized ultra-high pressure conditions. *Journal of Chromatography A* (2009) 1216, 4597-4605, Zhang, Y., Wang, X., Mukherjee, P. & Petersson, P.

References

54. Comparison study of porous, fused-core, and monolithic silica-based C18 HPLC columns for Celestoderm-V Ointment® analysis. *Journal of pharmaceutical and biomedical analysis* (2009) 50, 815-822, Zheng, J., Patel, D., Tang, Q., Markovich, R. & Rustum, A.
55. A high throughput dried DMSO Log D lipophilicity measurement based on 96-well shake-flask and atmospheric pressure photoionization mass spectrometry detection. *Journal of Chromatography A* (2010) 1217, 1950-1955, Alelyunas, Y., Pelosi-Kilby, L., Turcotte, P., Kary, M. & Spreen, R.
56. Halo Columns: New Generation Technology for High Speed Liquid Chromatography. *Journal of chromatographic science* (2010) 48, 386-394, Ali, I., Gaitonde, V. & Grahn, A.
57. How well does urinary lyso-Gb3 function as a biomarker in Fabry disease? *Clinica Chimica Acta* (2010) 411, 1906-1914, Auray-Blais, C., Ntwari, A., Clarke, J., Warnock, D., Oliveira, J., Young, S., Millington, D., Bichet, D., Sirrs, S. & West, M.
58. Physical characterization and evaluation of HPLC columns packed with superficially porous particles. *Journal of separation science* (2010) 33, 2547-2557, Baker, J., Vinci, J., Moore, A. & Colón, L.
59. The kinetic plot method applied to gradient chromatography: theoretical framework and experimental validation. *Journal of Chromatography A* (2010) 1217, 2787-2795, Broeckhoven, K., Cabooter, D., Lynen, F., Sandra, P. & Desmet, G.
60. Relationship between the Particle Size Distribution of Commercial Fully Porous and Superficially Porous High-Performance Liquid Chromatography Column Packings and their Chromatographic Performance. *Journal of Chromatography A* (2010) 1217, 7074-7081, Cabooter, D., Fanigliulo, A., Bellazzi, G. & Allieri, B.
61. Efficiency of the same neat silica column in hydrophilic interaction chromatography and per aqueous liquid chromatography. *Journal of Chromatography A* (2010) 1217, 683-688, Gritti, F., dos Santos Pereira, A., Sandra, P. & Guiochon, G.
62. Physical properties and structure of fine core-shell particles used as packing materials for chromatography: Relationships between particle characteristics and column performance. *Journal of Chromatography A* (2010) 1217, 3819-3843, Gritti, F., Leonardi, I., Abia, J. & Guiochon, G.
63. New trends in fast and high-resolution liquid chromatography: a critical comparison of existing approaches. *Analytical and Bioanalytical Chemistry* (2010) 397, 1069-1082, Guilleme, D., Ruta, J., Rudaz, S. & Veuthey, J.
64. Comparison of the mass transfer in totally porous and superficially porous stationary phases in liquid chromatography. *Analytical and Bioanalytical Chemistry* (2010) 397, 1307-1314, Kiss, I., Bacskaý, I., Kilr, F. & Felinger, A.
65. Extension of the C18 stationary phase knowledge by using the carotenoid test. *Journal of separation science* (2010) 33, 3097-3105, Lesellier, E.
66. Widespread but small-scale changes in the structural and dynamic properties of vaccinia virus poly (A) polymerase upon association with its processivity factor in solution. *Biochemistry* (2010) 49, 6247-6262, Li, C., Koter, M., Ye, X., Zhou, S., Chou, W., Luo, R. & Gershon, P.
67. Dissipation and Residues Detection of Diocetyl diethylenetriamine Cetate in Rice Plant and Environment by QuEChERS Method and Liquid Chromatography/Electrospray Tandem Mass Spectrometry. *Bulletin of environmental contamination and toxicology* (2010) 84, 596-601, Li, X., Jiang, Y., Shan, W. & Pan, C.
68. A quantitative investigation of fucosylated serum glycoproteins with application to esophageal adenocarcinoma. *Electrophoresis* (2010) 31, 1833-1841, Mann, B., Madera, M., Klouckova, I., Mechref, Y., Dobrolecki, L., Hickey, R., Hammoud, Z. & Novotny, M.
69. Instrumental considerations for the effective operation of short, highly efficient fused-core columns. Investigation of performance at high flow rates and elevated temperatures. *Journal of Chromatography A* (2010) 1217, 4561-4567, McCalley, D.
70. The use of charged aerosol detection with HPLC for the measurement of lipids. *Lipidomics, Methods in Molecular Biology* (2010) 579, 469-482, Plante, M., Bailey, B. & Acworth, I.
71. The Impact of Sampling Time on Peak Capacity and Analysis Speed in On-line Comprehensive Two-Dimensional Liquid Chromatography. *Journal of Chromatography A* (2010) 1217, 5700-5709, Potts, L., Stoll, D., Li, X. & Carr, P.
72. Subcellular distribution of okadaic acid in the digestive gland of *Mytilus galloprovincialis*: First evidences of lipoprotein binding to okadaic acid. *Toxicon* (2010) 55, 221-226, Rossignoli, A. & Blanco, J.
73. Wider Pore Superficially Porous Particles for Peptide Separations by HPLC. *Journal of chromatographic science* (2010) 48, 566-571, Schuster, S., Wagner, B., Boyes, B. & Kirkland, J.
74. Optimization and validation of a high performance liquid chromatography method for rapid determination of sinifloxacin, a novel fluoroquinolone in rat plasma using a fused-core C18-silica column. *Journal of pharmaceutical and biomedical analysis* (2010) 51, 889-893, Wang, S., Wen, J., Cui, L., Zhang, X., Wei, H., Xie, R., Feng, B., Wu, Y. & Fan, G.

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Global Chromatography Co., Ltd.